

HSD3B7 Antibody



PACO06208

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:20000, WB:1:500-1:2000

Protein Background:**Gene ID:**

HSD3B7

Uniprot

Q9H2F3

Synonyms:

HSD3B7; 3 beta-hydroxysteroid dehydrogenase type 7; 3 beta-hydroxysteroid dehydrogenase type VII; 3-beta-HSD VII; 3-beta-hydroxy-Delta(5)-C27 steroid oxidoreductase; C(27) 3-beta-HSD; Cholest-5-ene-3-beta; 7-alpha-diol 3-beta-dehydrogenase


Immunogen:

Synthesized peptide derived from the Internal region of human 3 β -HSD7.


Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of MCF7 HELA cells using 3 beta -HSD7 Polyclonal Antibody.



Western Blot analysis of COLO205 cells using 3 beta -HSD7 Polyclonal Antibody.