

ABCA8 Antibody



PACO06215

Product Information

Size:

50ug

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:40000, WB:1:500-1:2000

Protein Background:

Gene ID:

ABCA8

Uniprot

O94911

Synonyms:

ABCA8; KIAA0822; ATP-binding cassette sub-family A member 8


Immunogen:

Synthesized peptide derived from the Internal region of human ABCA8.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of HuvEc PC12 cells using ABCA8 Polyclonal Antibody.



Western Blot analysis of HuvEc cells using ABCA8 Polyclonal Antibody.