

ABCC13 Antibody



PACO06216

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

ELISA:1:40000, WB:1:500-1:2000,
IHC:1:100-1:300

Protein Background:**Gene ID:**

ABCC13

Uniprot

Q9NSE7

Synonyms:

ABCC13; C21orf73; PRED6; Putative ATP-binding cassette sub-family C member 13


Immunogen:

Synthesized peptide derived from the Internal region of human ABCC13.


Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of HELA cells using ABCC13 Polyclonal Antibody.



Western Blot analysis of HT29 cells using ABCC13 Polyclonal Antibody.