

AKR1A1 Antibody



PACO06242

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:40000, WB:1:500-1:2000

Protein Background:

Gene ID:

AKR1A1

Uniprot

P14550

Synonyms:

AKR1A1; ALDR1; ALR; Alcohol dehydrogenase [NADP(+)]; Aldehyde reductase; Aldo-keto reductase family 1 member A1


Immunogen:

Synthesized peptide derived from the C-terminal region of human AKR1A1.


Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of HELA HuvEc MCF7 cells using AKR1A1 Polyclonal Antibody.



Western Blot analysis of NIH-3T3 cells using AKR1A1 Polyclonal Antibody.