

APLP2 Antibody



PACO06257

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

ELISA:1:40000, WB:1:500-1:2000,
IHC:1:100-1:300

Protein Background:**Gene ID:**

APLP2

Uniprot

Q06481

Synonyms:

APLP2; APPL2; Amyloid-like protein 2; APLP-2; APPH; Amyloid protein homolog; CDEI box-binding protein; CDEBP


Immunogen:

Synthesized peptide derived from the Internal region of human APLP-2.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of 293 NIH-3T3 cells using APLP-2 Polyclonal Antibody.



Western Blot analysis of 3T3 cells using APLP-2 Polyclonal Antibody.