

ARHGAP18 Antibody



PACO06261

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

ELISA:1:40000, WB:1:500-1:2000,
IHC:1:100-1:300

Protein Background:

Gene ID:

ARHGAP18

Uniprot

Q8N392

Synonyms:

ARHGAP18; Rho GTPase-activating protein 18; MacGAP; Rho-type GTPase-activating protein 18


Immunogen:

Synthesized peptide derived from the C-terminal region of human ARHGAP18.


Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of NIH-3T3 cells using ARHGAP18 Polyclonal Antibody.



Western Blot analysis of Jurkat cells using ARHGAP18 Polyclonal Antibody.