

BCLAF1 Antibody



PACO06301

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:10000, WB:1:500-1:2000

Protein Background:

Gene ID:

BCLAF1

Uniprot

Q9NYF8

Synonyms:

BCLAF1; BTF; KIAA0164; Bcl-2-associated transcription factor 1; Btf


Immunogen:

Synthesized peptide derived from the Internal region of human BTF.


Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of NIH-3T3 cells using BTF Polyclonal Antibody.



Western Blot analysis of A549 cells using BTF Polyclonal Antibody.