

CDH16 Antibody



PACO06310

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC, IF

Recommended dilutions:

ELISA:1:10000, WB:1:500-1:2000,
IHC:1:100-1:300, IF:1:200-1:1000

Protein Background:**Gene ID:**

CDH16

Uniprot

O75309

Synonyms:

CDH16; Cadherin-16; Kidney-specific cadherin; Ksp-cadherin


Immunogen:

Synthesized peptide derived from the C-terminal region of human Cadherin-16.


Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of HuvEc cells using Cadherin-16 Polyclonal Antibody.



Western Blot analysis of A549 cells using Cadherin-16 Polyclonal Antibody.