

# TNFSF9 Antibody



PACO06328

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, IF

**Recommended dilutions:**

ELISA:1:20000, IF:1:200-1:1000

**Protein Background:**

**Gene ID:**

TNFSF9

**Uniprot**

P41273

**Synonyms:**

TNFSF9; Tumor necrosis factor ligand superfamily member 9; 4-1BB ligand; 4-1BBL

**Immunogen:**

Synthesized peptide derived from the Internal region of human CD137L.

**Storage:**

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

## Product Images

---

N/A

N/A