

STK17B Antibody



PACO06435

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

ELISA:1:40000, WB:1:500-1:2000,
IHC:1:100-1:300

Protein Background:**Gene ID:**

STK17B

Uniprot

O94768

Synonyms:

STK17B; DRAK2; Serine/threonine-protein kinase 17B; DAP kinase-related apoptosis-inducing protein kinase 2


Immunogen:

Synthesized peptide derived from the Internal region of human DRAK2.


Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of Jurkat K562 cells using DRAK2 Polyclonal Antibody.



Western Blot analysis of brain cells using DRAK2 Polyclonal Antibody.