

DNM1 Antibody



PACO06440

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

ELISA:1:40000, WB:1:500-1:2000,
IHC:1:100-1:300

Protein Background:

Gene ID:

DNM1

Uniprot

Q05193

Synonyms:

DNM1; DNM; Dynamin-1


Immunogen:

Synthesized peptide derived from human Dynamin I around the non-phosphorylation site of S774.


Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of NIH-3T3 cells using Dynamin I Polyclonal Antibody.



Western Blot analysis of 3T3 cells using Dynamin I Polyclonal Antibody.