

EEF2 Antibody



PACO06450

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC, IF

Recommended dilutions:

ELISA:1:5000, WB:1:500-1:2000, IHC:1:100-1:300, IF:1:200-1:1000

Protein Background:

Gene ID:

EEF2

Uniprot

P13639

Synonyms:

EEF2; EF2; Elongation factor 2; EF-2


Immunogen:

Synthesized peptide derived from human EF-2 around the non-phosphorylation site of T56.


Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of various cells using EF-2 Polyclonal Antibody.



Western Blot analysis of VEC cells using EF-2 Polyclonal Antibody.