

HNRNPL Antibody



PACO06558

Product Information

Size:

50ug

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC, IF

Recommended dilutions:

ELISA:1:20000, WB:1:500-1:2000,
IHC:1:100-1:300, IF:1:200-1:1000

Protein Background:**Gene ID:**

HNRNPL

Uniprot

P14866

Synonyms:

HNRNPL; HNRPL; P/OKcl.14; Heterogeneous nuclear ribonucleoprotein L; hnRNP L


Immunogen:

Synthesized peptide derived from the Internal region of human hnRNP L.


Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of HELA cells using hnRNP L Polyclonal Antibody.



Western Blot analysis of HuvEc cells using hnRNP L Polyclonal Antibody.