

IL1B Antibody



PACO06583

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

ELISA:1:5000, WB:1:500-1:2000, IHC:1:100-1:300

Protein Background:

Gene ID:

IL1B

Uniprot

P01584

Synonyms:

IL1B; IL1F2; Interleukin-1 beta; IL-1 beta; Catabolin


Immunogen:

Synthesized peptide derived from the Internal region of human IL-1 β .


Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of L929 MOUSE-BRAIN cells using IL-1 beta Polyclonal Antibody.



Western Blot analysis of HeLa cells using IL-1 beta Polyclonal Antibody.