

# IRAK2 Antibody



PACO06591

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:**

ELISA:1:40000, WB:1:500-1:2000

**Protein Background:**

**Gene ID:**

IRAK2

**Uniprot**

O43187

**Synonyms:**

IRAK2; Interleukin-1 receptor-associated kinase-like 2; IRAK-2

**Immunogen:**


Synthesized peptide derived from the Internal region of human IRAK-2.

**Storage:**


Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

## Product Images

---



Western Blot analysis of Jurkat HuvEc cells using IRAK-2 Polyclonal Antibody.



Western Blot analysis of 293 cells using IRAK-2 Polyclonal Antibody.