

MERTK Antibody



PACO06636

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:5000, WB:1:500-1:2000

Protein Background:

Gene ID:

MERTK

Uniprot

Q12866

Synonyms:

MERTK; MER; Tyrosine-protein kinase Mer; Proto-oncogene c-Mer; Receptor tyrosine kinase MerTK


Immunogen:

Synthesized peptide derived from the N-terminal region of human MerTK.


Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of various cells using MerTK Polyclonal Antibody.



Western Blot analysis of HEPG2 cells using MerTK Polyclonal Antibody.