

MNT Antibody



PACO06648

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:40000, WB:1:500-1:2000

Protein Background:

Gene ID:

MNT

Uniprot

Q99583

Synonyms:

MNT; BHLHD3; ROX; Max-binding protein MNT; Class D basic helix-loop-helix protein 3; bHLHd3; Myc antagonist MNT; Protein ROX


Immunogen:

Synthesized peptide derived from the Internal region of human MNT.


Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of Jurkat cells using MNT Polyclonal Antibody.



Western Blot analysis of A549 cells using MNT Polyclonal Antibody.