

ABCC10 Antibody



PACO06652

Product Information

Size:

50ug

Reactivity:

Human, Monkey

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:40000, WB:1:500-1:2000

Protein Background:

Gene ID:

ABCC10

Uniprot

Q5T3U5

Synonyms:

ABCC10; MRP7; SIMRP7; Multidrug resistance-associated protein 7; ATP-binding cassette sub-family C member 10


Immunogen:

Synthesized peptide derived from the Internal region of human MRP7.


Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of Jurkat A549 cells using MRP7 Polyclonal Antibody.



Western Blot analysis of COS7 cells using MRP7 Polyclonal Antibody.