

MRPL3 Antibody



PACO06654

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:5000, WB:1:500-1:2000

Protein Background:

Gene ID:

MRPL3

Uniprot

P09001

Synonyms:

MRPL3; MRL3; RPML3; 39S ribosomal protein L3; mitochondrial; L3mt; MRP-L3


Immunogen:

Synthesized peptide derived from the Internal region of human MRP-L3.


Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of Jurkat cells using MRP-L3 Polyclonal Antibody.



Western Blot analysis of Jurkat cells using MRP-L3 Polyclonal Antibody.