

OR10G4 Antibody



PACO06709

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IF

Recommended dilutions:

ELISA:1:10000, WB:1:500-1:2000, IF:1:200-1:1000

Protein Background:

Gene ID:

OR10G4

Uniprot

Q8NGN3

Synonyms:

OR10G4; Olfactory receptor 10G4; Olfactory receptor OR11-278


Immunogen:

Synthesized peptide derived from the C-terminal region of human Olfactory receptor 10G4.


Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of various cells using Olfactory receptor 10G4 Polyclonal Antibody.



Western Blot analysis of MCF7 cells using Olfactory receptor 10G4 Polyclonal Antibody.