

ORC1 Antibody



PACO06715

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:5000, WB:1:500-1:2000

Protein Background:

Gene ID:

ORC1

Uniprot

Q13415

Synonyms:

ORC1; ORC1L; PARC1; Origin recognition complex subunit 1; Replication control protein 1


Immunogen:

Synthesized peptide derived from the Internal region of human ORC1.


Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of NIH-3T3 HeLa 293 cells using ORC1 Polyclonal Antibody.



Western Blot analysis of 3T3 cells using ORC1 Polyclonal Antibody.