

SRCIN1 Antibody



PACO06718

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:20000, WB:1:500-1:2000

Protein Background:

Gene ID:

SRCIN1

Uniprot

Q9C0H9

Synonyms:

SRCIN1; KIAA1684; P140; SRC kinase signaling inhibitor 1; SNAP-25-interacting protein; SNIP; p130Cas-associated protein; p140Cap


Immunogen:

Synthesized peptide derived from the Internal region of human p140Cap.


Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of 293T HELA cells using p140Cap Polyclonal Antibody.



Western Blot analysis of HeLa cells using p140Cap Polyclonal Antibody.