

# RNF144B Antibody



PACO06729

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human, Mouse

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:**

ELISA:1:40000, WB:1:500-1:2000

**Protein Background:****Gene ID:**

RNF144B

**Uniprot**

Q7Z419

**Synonyms:**

RNF144B; IBRDC2; P53RFP; E3 ubiquitin-protein ligase RNF144B; IBR domain-containing protein 2; RING finger protein 144B; p53-inducible RING finger protein

**Immunogen:**


Synthesized peptide derived from the Internal region of human p53RFP.

**Storage:**


Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

## Product Images

---



Western Blot analysis of RAT-MUSCLE HELA 293T AD293 cells using p53RFP Polyclonal Antibody.



Western Blot analysis of HeLa cells using p53RFP Polyclonal Antibody.