

TAF6L Antibody



PACO06733

Product Information

Size:

50ug

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

ELISA:1:40000, WB:1:500-1:2000,
IHC:1:100-1:300

Protein Background:

Gene ID:

TAF6L

Uniprot

Q9Y6J9

Synonyms:

TAF6L; PAF65A; TAF6-like RNA polymerase II p300/CBP-associated factor-associated factor 65 kDa subunit 6L; PCAF-associated factor 65-alpha; PAF65-alpha


Immunogen:

Synthesized peptide derived from the N-terminal region of human PAF65 α .


Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of various cells using PAF65 alpha Polyclonal Antibody.



Western Blot analysis of 293 cells using PAF65alpha Polyclonal Antibody.