

EHHADH Antibody



PACO06737

Product Information

Size:

50ug

Reactivity:

Human, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

ELISA:1:10000, WB:1:500-1:2000,
IHC:1:100-1:300

Protein Background:

Gene ID:

EHHADH

Uniprot

Q08426

Synonyms:

EHHADH; ECHD; Peroxisomal bifunctional enzyme; PBE; PBEF


Immunogen:

Synthesized peptide derived from the Internal region of human PBEF.


Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of Jurkat K562 NIH-3T3 cells using PBFE Polyclonal Antibody.



Western Blot analysis of A549 cells using PBFE Polyclonal Antibody.