

MPRIP Antibody



PACO06802

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:5000, WB:1:500-1:2000

Protein Background:

Gene ID:

MPRIP

Uniprot

Q6WCQ1

Synonyms:

MPRIP; KIAA0864; MRIP; RHOIP3; Myosin phosphatase Rho-interacting protein; M-RIP; Rho-interacting protein 3; RIP3; p116Rip


Immunogen:

Synthesized peptide derived from the Internal region of human M-RIP.


Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of K562 Jurkat RAT-MUSCLE cells using M-RIP Polyclonal Antibody.



Western Blot analysis of 293 cells using M-RIP Polyclonal Antibody.