

ST5 Antibody



PACO06823

Product Information

Size:

50ug

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:20000, WB:1:500-1:2000

Protein Background:

Gene ID:

ST5

Uniprot

P78524

Synonyms:

ST5; DENND2B; HTS1; Suppression of tumorigenicity 5 protein; DENN domain-containing protein 2B; HeLa tumor suppression 1


Immunogen:

Synthesized peptide derived from the Internal region of human ST5.


Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of K562 cells using ST5 Polyclonal Antibody.



Western Blot analysis of COLO205 cells using ST5 Polyclonal Antibody.