

STK33 Antibody



PACO06827

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

ELISA:1:40000, WB:1:500-1:2000,
IHC:1:100-1:300

Protein Background:

Gene ID:

STK33

Uniprot

Q9BYT3

Synonyms:

STK33; Serine/threonine-protein kinase 33


Immunogen:

Synthesized peptide derived from the N-terminal region of human STK33.


Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of HuvEc Jurkat cells using STK33 Polyclonal Antibody.



Western Blot analysis of 293 cells using STK33 Polyclonal Antibody.