

VIM Antibody



PACO06855

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:10000, WB:1:500-1:2000

Protein Background:

Gene ID:

VIM

Uniprot

P08670

Synonyms:

VIM; Vimentin


Immunogen:

Synthesized peptide derived from human Vimentin around the non-phosphorylation site of Y61.


Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of 293 A549 Jurkat cells using Vimentin Polyclonal Antibody.



Western Blot analysis of Jurkat cells using Vimentin Polyclonal Antibody.