

CD53 Antibody



PACO06889

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:10000, WB:1:500-1:2000

Protein Background:

Gene ID:

CD53

Uniprot

P19397

Synonyms:

CD53; MOX44; TSPAN25; Leukocyte surface antigen CD53; Cell surface glycoprotein CD53; Tetraspanin-25; Tspan-25; CD53


Immunogen:

Synthesized peptide derived from the Internal region of human CD53.


Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of KB cells using CD53 Polyclonal Antibody.



Western Blot analysis of KB cells using CD53 Polyclonal Antibody.