

GPR119 Antibody



PACO06944

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IF

Recommended dilutions:

ELISA:1:10000, WB:1:500-1:2000, IF:1:200-1:1000

Protein Background:

Gene ID:

GPR119

Uniprot

Q8TDV5

Synonyms:

GPR119; Glucose-dependent insulinotropic receptor; G-protein coupled receptor 119


Immunogen:

Synthesized peptide derived from the Internal region of human GPR119.


Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western blot analysis of SH-SY5Y mouse-kidney lysis using GPR119 antibody.



Western Blot analysis of COLO cells using GPR119 Polyclonal Antibody.