

FAM117B Antibody



PACO06946

Product Information

Size:

50ug

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:5000, WB:1:500-1:2000

Protein Background:

Gene ID:

FAM117B

Uniprot

Q6P1L5

Synonyms:

FAM117B; ALS2CR13; Protein FAM117B; Amyotrophic lateral sclerosis 2 chromosomal region candidate gene 13 protein

Immunogen:

Synthesized peptide derived from the C-terminal region of human ALS2CR13.


Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western Blot analysis of HELA cells using ALS2CR13 Polyclonal Antibody.



Western Blot analysis of HELA cells using ALS2CR13 Polyclonal Antibody.