

ZBTB45 Antibody



PACO06996

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

WB, IHC

Recommended dilutions:

WB:1:500-1:1000, IHC:1:200-1:500

Protein Background:

Gene ID:

ZBTB45

Uniprot

Q96K62

Synonyms:

Zinc finger and BTB domain-containing protein 45


Immunogen:

Recombinant Protein.


Storage:

PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol.

Product Images



Western blot analysis of 1) Hela, 2) Jurkat, 3) Mouse Brain, 4) Mouse Kidney, 5) Rat Brain using ZBTB45 Polyclonal Antibody.



Immunohistochemical analysis of paraffin-embedded Mouse Kidney Tissue using ZBTB45 Polyclonal Antibody.