

MAPK6 Antibody



PACO07004

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

WB

Recommended dilutions:

WB:1:500-1:1000, IHC:1:200-1:500

Protein Background:

Gene ID:

MAPK6

Uniprot

Q16659

Synonyms:

Mitogen-activated protein kinase 6


Immunogen:

Recombinant Protein.


Storage:

PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol.

Product Images



Western blot analysis of 1) HeLa, 2) 293T, 3) 3T3, 4) Rat Brain using ERK 3 Polyclonal Antibody.



Immunohistochemical analysis of paraffin-embedded human Breast carcinoma using ERK 3 Polyclonal Antibody.