

# MAPK7 Antibody



PACO07005

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

WB, IHC

**Recommended dilutions:**

WB:1:1000-1:2000, IHC:1:200-1:500

**Protein Background:**

**Gene ID:**

MAPK7

**Uniprot**

Q13164

**Synonyms:**

Mitogen-activated protein kinase 7

**Immunogen:**


Recombinant Protein.

**Storage:**


PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol.

## Product Images

---



Western blot analysis of 1) HeLa, 2) 293T, 3) 3T3, 4) Mouse Skeletal Muscle, 5) Rat Kidney, 6) Rat Skeletal Muscle using ERK 5 Polyclonal Antibody.



Immunohistochemical analysis of paraffin-embedded human Breast carcinoma using ERK 5 Polyclonal Antibody.