

MAPK7 Antibody



PACO07006

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

WB, IHC

Recommended dilutions:

WB:1:1000-1:2000, IHC:1:200-1:500

Protein Background:

Gene ID:

MAPK7

Uniprot

Q13164

Synonyms:

Mitogen-activated protein kinase 7


Immunogen:

Recombinant Protein.


Storage:

PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol.

Product Images



Western blot analysis of 1) HeLa, 2) 293T, 3) Mouse Skeletal Muscle, 4) Rat Kidney, 5) Rat Skeletal Muscle using ERK 5 Polyclonal Antibody.



Immunohistochemical analysis of paraffin-embedded human Breast carcinoma using ERK 5 Polyclonal Antibody.