

Histone H2A. X Antibody



PACO07014

Product Information

Size:

50ug

Reactivity:

Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

WB, IHC

Recommended dilutions:

WB:1:1000-1:2000, IHC:1:200-1:500

Protein Background:

Gene ID:

H2AFX

Uniprot

P16104

Synonyms:

Histone H2AX


Immunogen:

Synthetic peptide.


Storage:

PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol.

Product Images



Western blot analysis of 1) 3T3, 2) Rat Brain Tissue using Histone H2A. X Polyclonal Antibody.



Immunohistochemical analysis of paraffin-embedded Mouse Colon Tissue using Histone H2A. X Polyclonal Antibody.