

CACNG3 Antibody



PACO07022

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

WB, IHC

Recommended dilutions:

WB:1:1000-1:2000, IHC:1:200-1:500

Protein Background:

Gene ID:

CACNG3

Uniprot

O60359

Synonyms:

Voltage-dependent calcium channel gamma-3 subunit


Immunogen:

Synthetic peptide.


Storage:

PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol.

Product Images



Western blot analysis of 1) Human Brain Tissue, 2) Mouse Brain Tissue, 3) Rat Brain Tissue using CACNG3 Polyclonal Antibody.



Immunohistochemical analysis of paraffin-embedded Rat Brain Tissue using CACNG3 Polyclonal Antibody.