

GNL3L Antibody




PACO07093


Product Information

Size: 50ug	Protein Background:
Reactivity: Human	Gene ID: GNL3L
Source: Rabbit	Uniprot Q9NVN8
Isotype: IgG	Synonyms: GNL3L; Guanine nucleotide-binding protein-like 3-like protein
Applications: ELISA, WB, IHC	Immunogen: Synthesized peptide derived from the Internal region of human GNL3L.
Recommended dilutions: ELISA:1:5000, WB:1:500-1:2000, IHC:1:100-1:300	Storage: Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.


Product Images



Western Blot analysis of NIH-3T3 LOVO cells using GNL3L Polyclonal Antibody.



Western blot analysis of mouse-brain HELA KB SH-SY5Y 293T 3T3 Mouse-kidney lysis using GNL3L antibody.



Western Blot analysis of HT29 cells using GNL3L Polyclonal Antibody.