

AWAT1 Antibody



PACO07132

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:5000, WB:1:500-1:2000

Protein Background:

Gene ID:

AWAT1

Uniprot

Q58HT5

Synonyms:

AWAT1; DGA2; DGAT2L3; Acyl-CoA wax alcohol acyltransferase 1; Diacylglycerol O-acyltransferase 2-like protein 3; Diacylglycerol acyltransferase 2; Long-chain-alcohol O-fatty-acyltransferase 1


Immunogen:

Synthesized peptide derived from the C-terminal region of human DGAT2L3.


Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.


Product Images



Western Blot analysis of 293T NIH-3T3 COLO cells using DGAT2L3 Polyclonal Antibody.



Western blot analysis of mouse-kidney KB lysis using DGAT2L3 antibody.



Western Blot analysis of 293 cells using DGAT2L3 Polyclonal Antibody.