

NOS3 Antibody



PACO07153

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IF

Recommended dilutions:

ELISA:1:10000, WB:1:500-1:2000, IF:1:200-1:1000

Protein Background:

Gene ID:

NOS3

Uniprot

P29474

Synonyms:

NOS3; Nitric oxide synthase; endothelial; Constitutive NOS; cNOS; EC-NOS; Endothelial NOS; eNOS; NOS type III; NOSIII


Immunogen:

Synthesized peptide derived from human NOS3 around the non-phosphorylation site of S1179.


Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.


Product Images



Western Blot analysis of A549 cells using NOS3 Polyclonal Antibody.



Western blot analysis of mouse-brain lysis using NOS3 antibody.



Western Blot analysis of A549 cells using NOS3 Polyclonal Antibody.