

CACNA1B Antibody



PACO07168

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC-p:1:50-1:100

Protein Background:

Gene ID:

CACNA1B

Uniprot

Q00975

Synonyms:

Immunogen:

Synthetic peptide.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Immunohistochemical analysis of paraffin-embedded Mouse Brain Tissue using Cav2.2Rabbit pAb diluted at 1:200.



Immunohistochemical analysis of paraffin-embedded Rat Brain Tissue using Cav2.2Rabbit pAb diluted at 1:200.