

GALR2 Antibody




PACO07179

Product Information

Size: 50ug	Protein Background:
Reactivity: Human, Mouse, Rat	Gene ID: GALR2
Source: Rabbit	Uniprot O43603
Isotype: IgG	Synonyms: GALR2; GALNR2; Galanin receptor type 2; GAL2-R; GALR-2
Applications: ELISA, WB	Immunogen: Synthetic peptide.
Recommended dilutions: ELISA:1:2000-1:10000, WB:1:500-1:1000	Storage: Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western blot analysis of 1) Rat Brain Tissue, 2) Mouse Brain Tissue with GALR2 Rabbit pAb diluted at 1:2,000.