

TRMT11 Antibody



PACO07209

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:10000, WB:1:500-1:2000

Protein Background:**Gene ID:**

TRMT11

Uniprot

Q7Z4G4

Synonyms:

TRMT11; C6orf75; MDS024; tRNA; guanine(10)-N2)-methyltransferase homolog; tRNA guanosine-2'-O-methyltransferase TRM11 homolog


Immunogen:

Synthesized peptide derived from the C-terminal region of human TRMT11.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images



Western blot analysis of HepG2 JK Colo HeLa using TRMT11 antibody. Secondary antibody was diluted at 1:20000.