

# CYP26B1 Antibody



PACO07288

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:**

ELISA:1:10000, WB:1:500-1:10000

**Protein Background:**

**Gene ID:**

CYP26B1

**Uniprot**

Q9NR63

**Synonyms:**

CYP26B1 CYP26A2 P450RAI2

**Immunogen:**

Synthetic peptide from human protein at AA range: 391-440.

**Storage:**

PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol.

## Product Images

---



Western Blot analysis of mouse-brain cells using Antibody diluted at 800. Secondary antibody was diluted at 1:20000.