

# GAB1 Antibody



PACO07309

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, IHC

**Recommended dilutions:**

ELISA:1:2000-1:10000, IHC-p:1:50-1:200

**Protein Background:**

**Gene ID:**

GAB1

**Uniprot**

Q13480

**Synonyms:**

GRB2-associated-binding protein 1 (GRB2-associated binder 1) (Growth factor receptor bound protein 2-associated protein 1)

**Immunogen:**


Recombinant Protein of Gab1.

**Storage:**

PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol.

## Product Images

---



Immunohistochemical analysis of paraffin-embedded Human Breast Carcinoma Tissue using Gab1 Rabbit pAb diluted at 1:500.