

IL17C Antibody




PACO07387


Product Information

Size:	Protein Background:
50ug	
Reactivity:	Gene ID:
Human	IL17C
Source:	Uniprot
Rabbit	Q9P0M4
Isotype:	Synonyms:
IgG	Interleukin-17C (IL-17C) (Cytokine CX2)
Applications:	Immunogen:
ELISA, IHC	Synthetic peptide from human protein at AA range: 81-130.
Recommended dilutions:	Storage:
ELISA:1:10000-1:20000, IHC-p:1:50-1:300	PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol.


Product Images



Immunohistochemical analysis of paraffin-embedded human-colon, antibody was diluted at 1:200.



Immunohistochemical analysis of paraffin-embedded human-pancreas, antibody was diluted at 1:200.



Immunohistochemical analysis of paraffin-embedded human-colon, antibody was diluted at 1:200.