

NOG Antibody




PACO07393

Product Information

| | |
|--|---|
| Size: 50ug | Protein Background: |
| Reactivity: Human, Mouse, Rat | Gene ID: NOG |
| Source: Rabbit | Uniprot Q13253 |
| Isotype: IgG | Synonyms: Noggin |
| Applications: ELISA, IHC | Immunogen: Synthetic peptide from human protein at AA range: 21-70. |
| Recommended dilutions: ELISA:1:10000-1:20000, IHC-p:1:50-1:300 | Storage: PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol. |

Product Images



Immunohistochemical analysis of paraffin-embedded human-liver, antibody was diluted at 1:200.