

ATG4C Antibody



PACO07469

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:10000, IHC:1:50-1:200

Protein Background:**Gene ID:**

ATG4C

Uniprot

Q96DT6

Synonyms:

Cysteine protease ATG4C (EC 3.4.22) (AUT-like 3 cysteine endopeptidase) (Autophagin-3) (Autophagy-related cysteine endopeptidase 3) (Autophagy-related protein 4 homolog C)


Immunogen:

Recombinant Protein.


Storage:

PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol.

Product Images



Immunohistochemical analysis of paraffin-embedded Human Brain Tissue using ATG4c Rabbit pAb diluted at 1:200.



Immunohistochemical analysis of paraffin-embedded Human Skeletal Muscle Tissue using ATG4c Rabbit pAb diluted at 1:200.