

ABCB6 Antibody



PACO07519

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

ABCB6

Uniprot

Q9NP58

Synonyms:

ATP-binding cassette, sub-family B (MDR/TAP), member 6; ABCB6; ABC; ABC14; EST45597; FLJ22414; MTABC3; PRP; umat ;

Immunogen:

Human ABCB6.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A